

600 HCN RT z-axis Probe Default Parameters

- **Tuning for Ubiquitin**
 - ^1H ~ 55
 - ^{13}C ~ 65
 - ^{15}N ~ 65
- **Offsets**
 - tof ~ 212 (~4.77 ppm at 25 deg C)
 - dof2 (^{15}N) ~ xxxx (~119 ppm)
- **Pulse widths**
 - ^1H : pw = 9.2 at tpwr = 59
 - ^{13}C : pwC = 12.3 at pwClvl = 61
 - ^{15}N : pwN = 41 at pwNlvl = 60
 - ^2H : H2pwD = xxxx at H2pwDlvl = xxxx
- **Deuterium decoupling using channel 4 (assumes garp1 or waltz with pwD90 of 250 usec)**
 - H2dpwr3D = xxxx
 - H2dmf3D = xxxx
- **Amplifier compressions**
 - compH = 1.00
 - compC = 0.90
 - compN = xxx
- **Maximum Power Limits (set in ghn_co)**
 - Channel 1: 60 dB
 - Channel 2: 61 dB
 - Channel 3: 61 dB
 - Channel 4: 50 dB
- **BioPack Power Limits (set in ghn_co)**
 - BPpwrlimits = 0
- **Gradients**
 - gzlvl1 = 30000, gzlvl2 = 30400 (ratio 0.987)
- Lock level values for lineshape sample
 - Lock power = 20
 - Lock gain = 23
 - Lock level when shimmed well > 80
- **Gain**
 - Gain values should always be > 18. Values between 20 and 25 are typically good values.